



ALIF Cage



- PEEK OPTIMA®
- User friendly instrumentation
- Titanium Markers for accurate positioning
- Variety of sizes to suit all your requirements
- Robotic Finger Type teeth minimize risk of implant movement

Unique Marketing through _____
innovation, design and
_____ manufacture.

Enquiries: +27 41 453 4605; info@ortho-sol.com

AVAILABLE SIZES:

Cage Code	Posterior Height (mm)	Anterior Height (mm)	Width (mm)	Depth (mm)	Instrument code
ALIF0406D18W25	4	6	25	18	OS4/6/18/25
ALIF0507D18W25	5	7	25	18	OS5/7/18/25
ALIF0608D18W25	6	8	25	18	OS6/8/18/25
ALIF0709D18W25	7	9	25	18	OS7/9/18/25
ALIF0810D18W25	8	10	25	18	OS8/10/18/25
ALIF0911D18W25	9	11	25	18	OS9/11/18/25
ALIF1012D18W25	10	12	25	18	OS10/13/18/25
ALIF1113D18W25	11	13	25	18	OS11/13/18/25
ALIF0406D21W30	4	6	30	21	OS4/6/21/30
ALIF0507D21W30	5	7	30	21	OS5/7/21/30
ALIF0608D21W30	6	8	30	21	OS6/8/21/30
ALIF0709D21W30	7	9	30	21	OS7/9/21/30
ALIF0810D21W30	8	10	30	21	OS8/10/21/30
ALIF0911D21W30	9	11	30	21	OS9/11/21/30
ALIF1012D21W30	10	12	30	21	OS10/12/21/30
ALIF1113D21W30	11	13	30	21	OS11/13/21/30
ALIF1214D21W30	12	14	30	21	OS12/13/21/30
ALIF1315D21W30	13	15	30	21	OS13/15/21/30

INSTRUMENTS:

- Easy to use instruments are available, comprising the following:
- “T” handle to hold the spacer trial. (Code: OSASTH)
- “T” handle for cage insertion (Code: OSACH)
- Spacer trials for each of the sizes in the table above (Codes as per table)

CAGE INSERTION TECHNIQUE

- Insert the trial into the space where the disc has been removed.
- Release the retractors.
- Ensure that the trial fits correctly in the disc space.
- If the fit is correct, prepare the cage size that corresponds to the trial size and prepare it with bone (allograft, or autograft bone). See appendix A.
- Insert the cage into the disc space, ensuring that it is positioned centrally with reference to the lateral aspects.
- The depth of the cage is determined by the stopper on the inserter. (This depth is ideally about 1mm posterior to the anterior cortex.)

APPENDIX A:

- The cage can be filled with either autograft bone (e.g. from the iliac crest), allograft bone (e.g. lyophilized bone, freeze dried bone), or else bone substitute (e.g. tri-calcium phosphate, bi-calcium phosphate, etc).
- The choice of void filler used is the surgeon’s preference after the necessary benefits have been carefully considered.
- Void filler should ideally protrude fractionally from both the cephalad and the caudal planes of the cage to allow for the initial resorption that occurs