



TRI-FIT CEMENTED STEM



- Highly Polished CoCr stem
- Rotational, axial and lateral stability
- Three tapered stem for maximum cement stability
- Tapered neck for optimum ROM and reduced impingement

Unique Marketing through _____
innovation, design and
_____manufacture.

Enquiries: +27 41 453 4605; info@ortho-sol.com

THE SATURN TRI-FIT HIP HAS BEEN DESIGNED BY A NUMBER OF LEADING ORTHOPAEDIC SURGEONS AND HAS BEEN MANUFACTURED USING ADVANCED ENGINEERING PRINCIPLES AND STATE-OF-THE-ART MACHINERY. THIS TOGETHER WITH CUTTING EDGE MATERIALS ENSURE THAT THE HIP IS OF THE HIGHEST QUALITY AVAILABLE IN ORDER TO BENEFIT THE PATIENT. A VARIETY OF SIZES ARE AVAILABLE TO FIT MOST PATIENT'S ANATOMY.

Stem Size	Offset	Catalogue code
1	45mm	638011
2	45mm	638012
3	45mm	638013
4	45mm	638014
1	38mm	638001
2	38mm	638002
3	38mm	638003
4	38mm	638004
CDH	36mm	638009

INDICATIONS

The Tri-Fit Stem is meant for use in total hip arthroplasties where a cemented stem is required, *inter alia* for the following:

- Ankylosis
- Failed previous hip surgery
- Avascular necrosis of the femoral neck
- Severe osteoarthritis or rheumatoid arthritis
- Acute traumatic fracture of the neck of the femur

PLEASE NOTE:

Total hip arthroplasty is restricted to surgeons that have undergone adequate training in this type of surgery.

The tri-Fit system should not be used with components other than those supplied by Ortho Sol Development (Pty) Ltd, as it may result in failure of the system

A centraliser and canal occluder are supplied with the Tri-Fit stem